



MARSHALL SURGICO
CORPORATION

PRODUCTS CATALOG
Complete Range of Electrosurgical Instruments



Symbol of Quality



info@marshallsurgico.com

www.marshallsurgico.com





Coaptation Nadler

L= 10.2 cm

Tip 0.5 mm



10-010-05

Iris (Str.)

L= 10.2 cm

Tip 0.5 mm



10-011-05

Iris (crv.)

L= 10.8 cm

Tip 1.0 mm



10-012-05

Straight

L= 16.0 cm

Tip 1.0 mm



10-016-05

Straight

L= 18.0 cm

Tip 1.0 mm



10-018-05

Straight

L= 20.0 cm

Tip 1.0 mm



10-020-05

Straight

L= 22.0 cm

Tip 1.0 mm



10-022-05

European Bipolar Forceps

Straight (Ang.)

L= 16.0 cm

Tip 1.0 mm



10-026-05

Straight (Ang.)

L= 18.0 cm

Tip 1.0 mm



10-028-05

Straight (Ang.)

L= 20.0 cm

Tip 1.0 mm



10-030-05

04

Straight (Ang.)

L= 22.0 cm

Tip 1.0 mm



10-032-05

Bayonet

L= 16.0 cm

Tip 1.0 mm



10-036-05

Bayonet

L= 18.0 cm

Tip 1.0 mm



10-038-05

Bayonet

L= 20.0 cm

Tip 1.0 mm



10-040-05



Bayonet (Ang.Up)

L= 16.0 cm

Tip 1.0 mm



10-046-10

Bayonet (Ang.Up)

L= 18.0 cm

Tip 1.0 mm



10-048-10

Bayonet (Ang.Up)

L= 20.0 cm

Tip 1.0 mm



10-050-10

Bayonet (Ang.Down)

L= 16.0 cm

Tip 1.0 mm



10-056-10

Bayonet (Ang.Up)

L= 18.0 cm

Tip 1.0 mm



10-058-10

Bayonet (Ang.Up)

L= 20.0 cm

Tip 1.0 mm



10-060-10

McPherson (Str.)

L= 9.0 cm

Tip 0.5 mm



101-001-05

European Bipolar Forceps

McPherson (crv.)

L= 11.5 cm

Tip 0.5 mm



101-002-05

Jeweler # 5

L= 11.5 cm

Tip 0.5 mm



101-003-05

Jeweler # 5

L= 11.5 cm

Tip 0.5 mm



101-003-05

06

Jeweler # 7

L= 11.5 cm

Tip 0.5 mm



101-008-05

Adson

L= 12.1 cm

Tip 1.0 mm



101-006-05

Straight

L= 16.0 cm

Tip 1.0 mm



101-016-05

Straight

L= 18.0 cm

Tip 1.0 mm



101-018-05



Straight

L= 20.0 cm

Tip 1.0 mm



101-020-10

Straight

L= 22.0 cm

Tip 1.0 mm



101-022-10

Straight (Ang.)

L= 16.0 cm

Tip 1.0 mm



101-028-10

Straight (Ang.)

L= 18.0 cm

Tip 1.0 mm



101-028-10

Straight (Ang.)

L= 20.0 cm

Tip 1.0 mm



101-030-10

Bayonet

L= 16.0 cm

Tip 1.0 mm



101-036-10

Bayonet

L= 18.0 cm

Tip 1.0 mm



101-038-10

European Bipolar Forceps

Bayonet

L= 22.0 cm

Tip 1.0 mm



Bayonet (Ang.Up)

L= 16.0 cm

Tip 1.0 mm



Bayonet (Ang.Up)

L= 18.0 cm

Tip 1.0 mm



08

Bayonet (Ang.Up)

L= 20.0 cm

Tip 1.0 mm



Bayonet (Ang.Down)

L= 16.0 cm

Tip 1.0 mm



Bayonet (Ang.Up)

L= 18.0 cm

Tip 1.0 mm



Bayonet (Ang.Up)

L= 20.0 cm

Tip 1.0 mm





European Irrigation Bipolar Forceps

Straight

L= 16.0 cm

Tip 1.0 mm



101-018-10 IR

Straight

L= 20.0 cm

Tip 1.0 mm



101-020-10 IR

Bayonet

L= 16.0 cm

Tip 1.0 mm



101-036-10 IR

Bayonet

L= 20.0 cm

Tip 1.0 mm



101-040-10 IR



European Non-Stick Bipolar Forceps

McPherson (Str.)

L= 9.0 cm

Tip 0.5 mm



10-001-05 NS

McPherson (crv.)

L= 11.5 cm

Tip 0.5 mm




10-002-05 NS

Jeweler # 5

L= 11.5 cm

Tip 0.5 mm



10-003-05 NS

10

Adson

L= 12.1 cm

Tip 1.0 mm



10-006-10 NS

Jeweler # 4

L= 11.5 cm

Tip 0.5 mm



10-007-05 NS

Jeweler # 7

L= 11.5 cm

Tip 0.5 mm




10-008-05 NS

Coaptation

L= 10.8 cm

Tip 0.5 mm



10-009-05 NS



European Non-Stick Bipolar Forceps

Coaptation Nadler

L= 10.2 cm

Tip 0.5 mm



10-010-05 NS

Iris (Str.)

L= 10.8 cm

Tip 0.5 mm



10-011-10 NS

Iris (crv.)

L= 10.8 cm

Tip 1.0 mm



10-012-10 NS

Straight

L= 18.0 cm

Tip 1.0 mm




10-018-10 NS

Straight

L= 20.0 cm

Tip 1.0 mm



10-019-10 NS

Straight

L= 20.0 cm

Tip 1.0 mm



10-020-10 NS

Straight

L= 22.0 cm

Tip 1.0 mm



10-022-10 NS



European Non-Stick Bipolar Forceps

Straight (Ang.)

L= 16.0 cm

Tip 1.0 mm



10-026-10 NS

Straight (Ang.)

L= 18.5 cm

Tip 1.0 mm



10-028-10 NS

Straight (Ang.)

L= 20.0 cm

Tip 1.0 mm



10-030-10 NS

12

Straight (Ang.)

L= 22.0 cm

Tip 1.0 mm



10-032-10 NS

Bayonet

L= 16.0 cm

Tip 1.0 mm



10-036-10 NS

Bayonet

L= 18.0 cm

Tip 1.0 mm



10-038-10 NS

Bayonet

L= 20.0 cm

Tip 1.0 mm



10-040-10 NS



European Non-Stick Bipolar Forceps

Bayonet (Ang.Up)

L= 16.0 cm

Tip 1.0 mm



10-046-10 NS

Bayonet (Ang.Up)

L= 18.0 cm

Tip 1.0 mm



10-048-10 NS

Bayonet (Ang.Up)

L= 20.0 cm

Tip 1.0 mm



10-050-10 NS

13

Bayonet (Ang.Down)

L= 16.0 cm

Tip 1.0 mm



10-056-10 NS

Bayonet (Ang.Up)

L= 18.0 cm

Tip 1.0 mm



10-058-10 NS

Bayonet (Ang.Up)

L= 20.0 cm

Tip 1.0 mm



10-060-10 NS

USA 2 Pin Bipolar Forceps

McPherson (Str.)

L= 9.0 cm

Tip 0.5 mm



08-001-05

McPherson (Crv.)

L= 9.0 cm

Tip 0.5 mm



08-002-05

Coaptation Nadler

L= 10.8 cm

Tip 0.5 mm



08-003-05

14

Coaptation

L= 10.8 cm

Tip 0.5 mm



08-004-05

Jeweler # 5

L= 11.0 cm

Tip 0.5 mm



08-005-05

Jeweler # 4

L= 11.5 cm

Tip 0.5 mm



08-006-05

Tenzel








L= 11.5 cm

Tip 0.5 mm



08-007-05



Jeweler # 7		08-008-05
L= 11.5 cm		
Tip 0.5 mm		
Iris (Str.)		08-009-10
L= 10.8 cm		
Tip 0.5 mm		
Iris (Crv.)		08-010-10
L= 10.8 cm		
Tip 1.0 mm		
Adson		08-011-10
L= 12.1 cm		
Tip 1.0 mm		
Semkin		08-012-10
L= 14.0 cm		
Tip 1.0 mm		
Cushing		08-013-10
L= 16.0 cm		
Tip 1.0 mm		
Cushing		08-014-10
L= 11.5 cm		
Tip 1.0 mm		

USA 2 Pin Bipolar Forceps

Cushing (Gursh)

L= 18.0 cm

Tip 1.0 mm



08-018-10

Gerald (Str)

L= 17.8 cm

Tip 1.0 mm



08-019-10

Gerald (Crv)

L= 17.8 cm

Tip 1.0 mm



08-020-10

16

Gerald Bayonet

L= 17.8 cm

Tip 1.0 mm



08-021-10

Jansen Bayonet

L= 15.0 cm

Tip 1.0 mm



08-027-10

Jansen Bayonet

L= 18.0 cm

Tip 1.0 mm



08-028-10

Mells

L= 17.8 cm

Tip 1.0 mm



08-026-10



Scollive Green-Wood

L= 19.1 cm

Tip 1.5 mm



08-025-15

Yasargil

L= 19.8 cm

Tip 1.0 mm



08-022-10

Hardy Bayonet

L= 12.5 cm

Tip 1.0 mm



08-032-10

Hardy Bayonet (Ang)

L= 21.5 cm

Tip 1.0 mm



08-033-10



USA 2-Pin Non-Stick Bipolar Forceps

McPherson (Str)

L= 9.0 cm

Tip 0.5 mm



McPherson (Crv)

L= 9.0 cm

Tip 0.5 mm



Jewelery #5

L= 11.5 cm

Tip 0.5 mm



18

Jeweler # 4

L= 11.5 cm

Tip 0.5 mm



Adson

L= 12.1 cm

Tip 1.0 mm



Semkin

L= 14.0 cm

Tip 1.0 mm



Cushing

L= 16.0 cm

Tip 1.0 mm





USA 2-Pin Non-Stick Bipolar Forceps

Cushing

L= 18.0 cm

Tip 1.0 mm



08-014-10-NS

Gerald Bayonet

L= 17.8 cm

Tip 1.0 mm



08-021-10-NS

Jansen Bayonet

L= 18.0 cm

Tip 1.0 mm



08-028-10-NS

Scollive Green-Wood

L= 19.1 cm

Tip 1.5 mm



08-025-15-NS

Yesargil

L= 19.8 cm

Tip 1.0 mm



08-022-10-NS

Hardy Bayonet

L= 21.5 cm

Tip 1.0 mm



08-032-10-NS

USA 2 Pin Bipolar Forceps

Semkln
Irrigation



09-012-IR

Gerald Bayonet
Irrigation



09-020-IR

Jansen Bayonet
Irrigation



09-026-IR

Scollive Green-Wood
Irrigation



09-031-IR

Hardy Bayonet
Irrigation



09-037-IR



MARSHALL SURGICO
CORPORATION

SilverGlide Non-Stick Bipolar Forceps

McPherson (Str.)

L= 9.0 cm

Tip 0.5 mm



08-001-05FNS

McPherson (Crv.)

L= 9.0 cm

Tip 0.5 mm



08-002-05FNS

Jeweler # 5

L= 11.5 cm

Tip 0.5 mm



08-005-05FNS

Adson

L= 12.1 cm

Tip 1.0 mm



08-011-10FNS

Cushing (Str.)

L= 16.0 cm

Tip 1.0 mm



08-016-10FNS

Jansen Bayonet

L= 18.0 cm

Tip 1.0 mm



08-028-10FNS

Jansen Bayonet

L= 20.0 cm

Tip 1.0 mm



08-030-10FNS

Bayonet Full Non-Stick

L= 16.0 cm

Tip 1.0 mm



30-036-10FNS

21



SilverGlide Non-Stick Bipolar Forceps

Bayonet Full Non-Stick

L= 18.0 cm

Tip 1.0 mm



30-038-10FNS

Bayonet Full Non-Stick

L= 20.0 cm

Tip 1.0 mm



30-040-10FNS

Bayonet Full Non-Stick

L= 22.0 cm

Tip 1.0 mm



30-042-10FNS

Bayonet (Irrigation)

L= 20.0 cm

Tip 1.0 mm



10-040-10FNS IR

Bayonet (Irrigation)

L= 22.0 cm

Tip 1.0 mm



10-042-10FNS IR

Bayonet Full Non-Stick

L= 22.0 cm

Tip 1.0 mm



10-042-10FNS

Bayonet Full Non-Stick

L= 22.0 cm

Tip 0.5 mm



48-022-05 FNS

Bayonet Full Non-Stick

L= 24.0 cm

Tip 1.0 mm



58-024-10 FNS



SilverGlide Non-Stick Bipolar Forceps

Jeweler # 5

L= 11.5 cm

Tip 0.5 mm



10-005-05FNS

Jeweler # 5

L= 11.5 cm

Tip 0.5 mm



10-006-05FNS

Straight Full Non-Stick

L= 12.0 cm

Tip 1.0 mm



10-012-10FNS

Straight (Ang.)

L= 16.0 cm

Tip 1.0 mm



10-026-10FNS

23

Bayonet Full Non-Stick

L= 18.0 cm

Tip 1.0 mm



10-038-10FNS

Bayonet Full Non-Stick

L= 20.0 cm

Tip 1.5 mm



10-040-15FNS

Bayonet Full Non-Stick

L= 22.0 cm

Tip 1.5 mm



10-042-15FNS

Bayonet Full Non-Stick

L= 22.0 cm

Tip 1.0 mm



10-044-10FNS



Single Use USA 2-Pin Bipolar Forceps

Gerald Straight
Single Use



08-019S

Jansen Bayonet
Irrigation



08-026 IRS

Jansen Bayonet
Single Use



08-030S

24

Jeweler # 5
Single Use



08-005S

Adson Forceps
Single Use



08-011S

Cushing Forceps
Size: 16.0 cm
Tip: 1.0 mm



08-016S

Jansen Bayonet ang-up
Size: 20.0 cm
Tip: 1.0 mm



08-042S



MARSHALL SURGICO
CORPORATION

SilverGlide Non-Stick Bipolar Forceps

Jeweler # 5

L= 11.5 cm

Tip 0.5 mm



10-005-05FNS

Jeweler # 5

L= 11.5 cm

Tip 0.5 mm



10-006-05FNS

Straight Full Non-Stick

L= 12.0 cm

Tip 1.0 mm



10-012-10FNS

Straight (Ang.)

L= 16.0 cm

Tip 1.0 mm



10-026-10FNS

25

Bayonet Full Non-Stick

L= 18.0 cm

Tip 1.0 mm



10-038-10FNS

Bayonet Full Non-Stick

L= 20.0 cm

Tip 1.5 mm



10-040-15FNS

Bayonet Full Non-Stick

L= 22.0 cm

Tip 1.5 mm



10-042-15FNS

Bayonet Full Non-Stick

L= 22.0 cm

Tip 1.0 mm



10-044-10FNS



Single Use USA 2-Pin Bipolar Forceps

Jansen Bayonet

Size: 18.0 cm

Tip: 1.0 mm



08-028S

Jeweler # 4

Size: 11.5 cm

Tip: 0.5 mm



08-004S

26

Adson Forceps

Size: 12.1 cm

Tip: 1.0 mm



18-011S

Jansen Bayonet

Size: 18.0 cm

Tip: 1.0 mm



18-028S

Jansen Bayonet

Size: 20.0 cm

Tip: 1.0 mm



18-030S



MARSHALL SURGICO
CORPORATION

Single Use European Bipolar Forceps

Bayonet Forceps

Size: 18.0 cm

Tip: 1.0 mm



10-038S

Bayonet Forceps

Size: 20.0 cm

Tip: 1.0 mm



10-040S

Stragng Forceps

Size: 20.0 cm

Tip: 1.0 mm



10-020S

Straight Forceps (ang)

Size: 14.0 cm

Tip: 1.0 mm



10-024S

Straight Forceps

Size: 14.0 cm

Tip: 1.0 mm



10-014S

27



Single Use USA 2-Pin Bipolar Forceps

Jeweler # 5
Single Use



10-003S

Straight Forceps
Single Use



10-020S

28

Bayonet Forceps
Single Use



10-040S

Straight Forceps
Single Use



10-019S

Straight Forceps (ang)
Single Use



10-030S



Silicone Monopolar Cables 10" (3.0m)
SB-1501
Instrument Connector: Screw Type
Generator Connector: Standard Bovie 8.0mm Plug
Compatibility: Valleylab, Conmed, Bovie



Silicone Monopolar Cables 10" (3.0m)
SB-1502
Instrument Connector: 4.0mm Male Fitting
Generator Connector: Standard 4.0mm Banana Plug
Compatibility: Valleylab, Conmed & All ES Generators



Silicone Monopolar Endoscopic Cable 10" (3.0m)
SB-1503
Instrument Connector: 4.0mm Female Plug
Generator Connector: Standard Bovie 8.0mm Plug
Compatibility: Valleylab, Conmed, Bovie



Silicone Monopolar Endoscopic Cable 10" (3.0m)
SB-1504
Instrument Connector: 4.0mm Female Plug
Generator Connector: Standard 4.0mm Banana Plug
Compatibility: Valleylab, Conmed & All ES Generators



Silicone Diathermy Cable 10" (3.0m)
SB-1505
Instrument Connector: 4.8mm Male Fitting
Generator Connector: Standard Bovie 8.0mm Plug
Compatibility: Valleylab, Conmed, Bovie



Silicone Diathermy Cable 10" (3.0m)
SB-1506
Instrument Connector: 4.8mm Male Fitting
Generator Connector: Standard 4.0mm Banana Plug
Compatibility: Valleylab, Conmed & All ES Generators



Bipolar Forceps Cables



Silicone Bipolar Cables 10" (3.0m)
SB 2501
Instrument Connector: USA 2-Pin Plug
Generator Connector: 4.0mm Banana Plug
Compatibility: Valleylab, Conmed, Codman, Olympus and Mbst Bipolar Generator



Silicone Bipolar Cables 10" (3.0m)
SB 2502
Instrument Connector: European Flat Plug
Generator Connector: Martin Type
Compatibility: Aesculap, Martin, Wolf, Berchold



Silicone Bipolar Cables 10" (3.0m)
SB 2503
Instrument Connector: European Flat Plug
Generator Connector: 4.0mm Banana Plug
Compatibility: Valleylab, Conmed, Codman, Olympus & Mbst Bipolar Generator



Silicone Bipolar Cables 10" (3.0m)
SB 2504
Instrument Connector: USA2-Pin Plug
Generator Connector: Martin Type
Compatibility: Aesculap, Martin, Wolf, Berchold

30



Silicone Bipolar Cables 10" (3.0m)
SB 2606
Instrument Connector: USA2-Pin Plug
Generator Connector: Erbe Type
Compatibility: Erbe, Siemens, Downs



Silicone Bipolar Cables 10" (3.0m)
SB 2506
Instrument Connector: USA2-Pin Plug
Generator Connector: Pin Distance 29.0mm
Compatibility: Valleylab, EMS, Dolley

Silicone Bipolar Cables 10" (3.0m)
SB 2507
Instrument Connector: European Flat Plug
Generator Connector: Pin Distance 29.0mm
Compatibility: Valleylab, EMS, Dolley





**Tumor Warwick Scissor
(Straight)**
SB-05-001



**Tumor Warwick Scissor
(Curved)**
SB-05-002



**Surgical Scissor
(Straight)**
SB-05-011



**Surgical Scissor
(Curved)**
SB-05-012



Abbey Sub Needle
17.5cm
SB-05-015



**Willson Hey Forceps
(Straight)**
SB-05-020

**Richie Forceps
(Straight)**
SB-05-038



**Willson Hey Forceps
(Curved)**
SB-05-022



**Artery Forceps
(Straight)**
SB-05-028

**Artery Forceps
(Curved)**
SB-05-028



Diathermy Instruments

Adson



SB-06-001

Mcldona



SB-06-002

Gillette



SB-06-003

Dissecting Forceps



SB-06-004

32

Gerald (1x2)



SB-06-005

Adson Bayonet



SB-06-006

Adson Bayonet



Angle Up

SB-06-016

Remington Hob



SB-06-019



Artery Forceps BI-Clamp

SB-18-004-25



Artery Forceps Non Stick

SB-18-006-25




Artery Forceps

SB-17-004-25



Artery Forceps

SB-17-001-25



Silicone BI-Clamp Connecting cord

Bipolar Electrodes

34



Straight

SB-14-001

Angled

SB-14-002

Angled

SB-14-003

Angled 45°

SB-14-004

Angled 90°

SB-14-005



BI-Turbinate Probe

Size: 17.5cm

SB-14-028



Electrodes 2.4mm



Blade Electrode
Disposable: 03-001



Ball Electrode (3.0 mm)
Disposable: 03-007



Circle Loop Electrode 4.0 mm
Reusable: 04-001



Ball Electrode (crv.) (3.0 mm)
Disposable: 03-008



Blade Electrode (crv.)
Disposable: 03-002



Circle Loop Electrode 4.0 mm
Reusable: 04-001



Needle Electrode
Disposable: 03-003



Ball Electrode (4.0 mm)
Disposable: 03-009



Circle Loop Electrode 8.0 mm
Reusable: 04-003



Needle Electrode (crv.)
Disposable: 03-004



Ball Electrode (crv.) (4.0 mm)
Disposable: 03-010



Diamond Electrode 8.0 mm
Reusable: 04-002



Ball Electrode (2.0 mm)
Disposable: 03-005



Ball Electrode (5.0 mm)
Disposable: 04-011




Tip Cleaner 03-35



Ball Electrode (crv.) (2.0 mm)
Disposable: 03-006



Ball Electrode (crv.) (5.0 mm)
Disposable: 03-012



Holder 03-38



Blade Electrode (Shaft 13.0 cm)
Disposable: 03-013



Blade Electrode (Shaft 13.0 cm)
Disposable: 03-015



Blade Electrode (Shaft 13.0 cm)
Disposable: 03-023



Blade Electrode (Shaft 13.0 cm)
Disposable: 03-017



Coagulation Electrode
Disposable: 03-025

Conization Electrodes



Conization Electrode
Disposable #03-028



Conization Electrode
Disposable #03-027



Conization Electrode
Disposable #03-029



Conization Electrode
Disposable #03-030



Conization Electrode
Disposable #03-030



Extender Electrode
03-032 Size: 7.0 cm



Extender Electrode
#03-033 Size:13.0 cm



Extender Electrode **#03-034 Size:18.0 cm**

36



Blade Electrode
Reusable 20-001



Ball Electrode
Reusable 20-007



Blade Electrode (Crv)
Reusable 20-002



Ball Electrode (Crv) 3mm
Reusable 20-008



Needle Electrode
Reusable 20-003



Ball Electrode (4mm)
Reusable 20-009



Needle Electrode (Crv)
Reusable 20-004



Ball Electrode (Crv) (4mm)
Reusable 20-010



Ball Electrode (2mm)
Reusable 20-005



Ball Electrode (6mm)
Reusable 20-011



Ball Electrode (Crv) (2mm)
Reusable 20-006



Ball Electrode (Crv) (5mm)
Reusable 20-012



Blade Electrode (Shaft 13.0cm)
Reusable 20-018



Needle Electrode (Shaft 13.0cm)
Reusable 20-019



Ball Electrode 3mm (Shaft 13.0 cm)
Reusable 20-020



Ball Electrode 4mm (Shaft 13.0cm)
Reusable 20-021



Ball Electrode 5mm (Shaft 13.0cm)
Reusable 20-022



Round Loop Electrode
(7.0cm) 4.0mm
Reusable 20-013



Round Loop Electrode
(7.0cm) 6.0mm
Reusable 20-014



Round Loop Electrode
(7.0cm) 8.0mm
Reusable 20-015



Round Loop Electrode
(7.0cm) 10.0mm
Reusable 20-016



Plate Electrode Small
Reusable 201-018



Plate Electrode Large
Reusable 20-017

Patient Return Plate

**Patient Return Plate Reusable
5m Cable #25-001**



38



**Patient Return Plate Reusable
(5m Cable) 25-002**



**Patient Return Cable
25-006-a (5m cable)
25-006-b (3m Cable)**



**Patient Return Cable
25-007-a (5m cable)
25-007-b (3m Cable)**



**S.S Patient Return Plate
(16x12cm)25-003**



WITHOUT SUCTION



15-001 L-Hook Electrode (36cm)



15-002 J-Hook Electrode (36cm)



15-003 Needle-Hook Electrode (36cm)



15-004 Spatula-Hook Electrode (36cm)



15-005 Ball-Hook Electrode (36cm)

WITH SUCTION



15-006 L-Hook Electrode (Irrigation) (36cm)



15-007 J-Hook Electrode (Irrigation) (36cm)



15-008 Needle-Hook Electrode (Irrigation) (36cm)



15-009 Spatula-Hook Electrode (Irrigation) (36cm)



15-010 Ball-Hook Electrode (Irrigation) (36cm)

Laparoscopic Electrodes





MARSHALL SURGICO
CORPORATION

Electrosurgical Instruments for Gynecology



Graves Vaginal Speculum
with smoke evacuation tube
Small: 21-001
Medium: 21-002
Large: 21-003



VU-Max Speculum
with smoke evacuation tube
Small: 21-011
Medium: 21-012
Large: 21-013



Pederson Vaginal Speculum
with smoke evacuation tube
Small: 21-041
Medium: 21-042
Large: 21-043

40



VU-More (Graves) Speculum
with smoke evacuation tube
Small: 21-031
Medium: 21-032
Large: 21-033

Cusco Speculum
with smoke evacuation tube
Small: 21-021
Medium: 21-022
Large: 21-023



Collin Speculum
with smoke evacuation tube
Small: 21-051
Medium: 21-052
Large: 21-053





Electrosurgical Instruments for Gynecology



Welman Speculum with Suction (Left Open)
Small: 21-061
Medium: 21-062
Large: 21-063



Welman Speculum with Suction (Right Open)
Small: 21-071
Medium: 21-072
Large: 21-073



Dressing Forceps
23-001 (21.0 cm)
23-001 (25.0 cm)



Tissue Forceps
23-003 (21.0 cm)
23-004 (25.0 cm)



Iris Hook
23-005 (25.0 cm)



Emmet Hook
23-006 (25.0 cm)



Straight Hook
23-007 (25.0 cm)



Two-Prong Hook
23-008 (25.0 cm)



Electrosurgical Instruments for Gynecology



42

**Lateral Vaginal
Retractor - #22-001**



**Kogan Endocervical
Speculum - #22-003**



**Lateral Vaginal Retractor
(wide-Open) #22-006**



**Rigby Retractor
22-007**



**Rigby Retractor
22-005**



3mm 130-133



4mm 130-134



3mm 130-133



MARSHALL SURGICO
CORPORATION

Electrosurgical Instruments for Gynecology



**Kogan Endocervical
Speculum**

SB-#25-003



**Lateral Vaginal
Retractor**

SB-#25-001



**Graves Speculum
with smoking evacuation tube**

**Small: 21-001
Medium: 21-002
Large: 21-003**



**Pederson Graves Speculum
with smoking evacuation tube**

**Small: 21-001
Medium: 21-002
Large: 21-003**

**Collin Speculum
with smoking evacuation tube**

**Small: 24-051
Medium: 24-052
Large: 24-053**



43



MARSHALL SURGICO
CORPORATION

Daska Road ,Capricorn Street,
Daska 51010, Punjab, Pakistan

Cell: +92 3445577007

Cell: +92 3337559333

Email: info@marshallsurgico.com

Web: www.marshallsurgico.com